


HANAROSTENT®
Biliary (NNN)



Indication :

- Palliative treatment of bile duct stricture caused by malignant tumors

Features :

- Flare on the both ends for Anti-Migration
- Uncovered stent
- Radiopaque markers for fluoroscopic guidance

Endoscopic Application

	Stent (mm)			Delivery device	
	Diameter	Usable length	Total length*	Length(mm)	Diameter(mm/fr)
SHS-08-040-180	11.5-8-11.5	26	40	1800	2.33/7
SHS-08-120-180		106	120	1800	2.33/7
SHS-10-040-180	13.5-10-13.5	26	40	1800	2.33/7
SHS-10-120-180		106	120	1800	2.33/7

*Increment of 10mm

	(NNN) Non covered	<u>XXX-00-000-000</u>	<u>XXX 00 000-X000</u>
	(CCC) Fully covered	① ② ③ ④	① ② ③ ④
	(NCN) Partially covered	① Product Name	③ Stent Length
		② Stent Diameter	④ Delivery Device Length