



Sphincterotome Model	PK PPK	Standard Precut
Catheter Diameter	5 7 7/5 8/6	5 F 7 F Reduced From 7F to 5F Reduced From 8F to 6F
Catheter Type	2L 3L	Double-Lumen Triple-Lumen
Cutting Wire Length	20 25 30	20 mm 25 mm 30 mm
Catheter Length	12 18 21 23 25 31	120 cm 180 cm 215 cm 230 cm 250 cm 310 cm
Guidewire Diameter	25 35	0,025" 0,035"
Cutting Wire Type	R W	Rope Wire