



STERILE EO

Coding

KYK 5 20 12 SU 3
 1 2 3 4 5 6

- 1- **KYK:** Forceps Model
- 2- **5:** Catheter Diameter
- 3- **20:** Forceps Opening
- 4- **12:** Catheter Length
- 5- **SU:** Handle Feature
- 6- **3:** Forceps Type

Forceps Model	KYK HYK	Hook-end Ring-end
Catheter Diameter	5 7	5 F (1,8 mm) 7 F (2,3 mm)
Forceps Opening	20 25 30	20 mm 25 mm 30 mm
Catheter Length	12 18 21 23 25 31	120cm 180 cm 215 cm 230 cm 250 cm 310 cm
Handle Feature	SU SY RD RY	Fixed End Fixed End Flushing Channel Rotatable Rotatable Flushing Channel
Forceps Type	3 4	Tripod Tetrapod