


**HANAROSTENT®**

**Duodenum/Pylorus Lasso (NNN)**



**Indication :**

- Palliative treatment of pyloric or duodenal stricture caused by malignant tumors

**Features :**

- Flare heads on the both ends for Anti-Migration
- Double layered wire structure and high flexibility
- Uncovered stent and lasso on the proximal end

**Fluoroscopic Application**

	Stent (mm)			Delivery device	
	Diameter	Usable length	Total length*	Length(mm)	Diameter(mm/fr)
DNZL-20-060-180	25-20-25	30	60	1800	3.4/10.2
DNZL-20-170-180		140	170	1800	3.4/10.2
DNZL-22-060-180	26-22-26	30	60	1800	3.4/10.2
DNZL-22-170-180		140	170	1800	3.4/10.2

\*Increment of 10mm

**Fluoroscopic Application**


	Stent (mm)			Delivery device	
	Diameter	Usable length	Total length*	Length(mm)	Diameter(mm/fr)
DNZL-20-060-120	25-20-25	30	60	1200	3.4/10.2
DNZL-20-170-120		140	170	1200	3.4/10.2
DNZL-22-060-120	26-22-26	30	60	1200	3.4/10.2
DNZL-22-170-120		140	170	1200	3.4/10.2

\*Increment of 10mm

	(NNN) Non covered	<u>XXX-00-000-000</u>	<u>XXX 00 000-X000</u>
	(CCC) Fully covered	① ② ③ ④	① ② ③ ④
	(NCN) Partially covered	① Product Name	③ Stent Length
		② Stent Diameter	④ Delivery Device Length

**HANAROSTENT®**

**Duodenum/Pylorus Lasso (NNN)**



**Indication :**

- Palliative treatment of pyloric or duodenal stricture caused by malignant tumors

**Features :**

- Flare heads on the both ends for Anti-Migration
- Double layered wire structure and high flexibility
- Uncovered stent and lasso on the proximal end

**Endoscopic Application**

	Stent (mm)			Delivery device	
	Diameter	Usable length	Total length*	Length(mm)	Diameter(mm/fr)
DNZL-20-060-230	25-20-25	30	60	2300	3.4/10.2
DNZL-20-170-230		140	170	2300	3.4/10.2
DNZL-22-060-230	26-22-26	30	60	2300	3.4/10.2
DNZL-22-170-230		140	170	2300	3.4/10.2

\*Increment of 10mm

	(NNN) Non covered	<u>XXX-00-000-000</u>	<u>XXX 00 000-X000</u>
	(CCC) Fully covered	① ② ③ ④	① ② ③ ④
	(NCN) Partially covered	① Product Name	③ Stent Length
		② Stent Diameter	④ Delivery Device Length