


# Niti-S Biliary Covered Stent [Bumpy]

For benign and/or malignant biliary strictures and benign pancreatic ductal strictures



## Features

- **Easy removal:** Fully PTFE (body portion) and silicone (both flared ends) covered design along with a removal string at the proximal end of the stent for easy removal
- **Anti-migration:** Both flared ends prevent the risk of migration
- **Radiopaque marker:** 3 (three) at both ends & 2 (two) in the middle

# Ordering Information

Percutaneous				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
TK0604CW	6	4	8.5	50
TK0605CW		5		
TK0606CW		6		
TK0607CW		7		
TK0608CW		8		
TK0609CW		9		
TK0610CW		10		
TK0612CW		12		
TK0804CW	8	4	8.5	50
TK0805CW		5		
TK0806CW		6		
TK0807CW		7		
TK0808CW		8		
TK0809CW		9		
TK0810CW		10		
TK0812CW		12		
TK1004CW	10	4	8.5	50
TK1005CW		5		
TK1006CW		6		
TK1007CW		7		
TK1008CW		8		
TK1009CW		9		
TK1010CW		10		
TK1012CW		12		

Percutaneous				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
TK0804CWS	8	4	8.5	50
TK0805CWS		5		
TK0806CWS		6		
TK0807CWS		7		
TK0808CWS		8		
TK0809CWS		9		
TK0810CWS		10		
TK0812CWS		12		
TK1004CWS	10	4	8.5	50
TK1005CWS		5		
TK1006CWS		6		
TK1007CWS		7		
TK1008CWS		8		
TK1009CWS		9		
TK1010CWS		10		
TK1012CWS		12		

# Ordering Information

Endoscopic				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
BK0604CW	6	4	8.5	180
BK0605CW		5		
BK0606CW		6		
BK0607CW		7		
BK0608CW		8		
BK0609CW		9		
BK0610CW		10		
BK0612CW		12		
BK0804CW		8		
BK0805CW	5			
BK0806CW	6			
BK0807CW	7			
BK0808CW	8			
BK0809CW	9			
BK0810CW	10			
BK0812CW	12			
BK1004CW	10		4	8.5
BK1005CW		5		
BK1006CW		6		
BK1007CW		7		
BK1008CW		8		
BK1009CW		9		
BK1010CW		10		
BK1012CW		12		

Endoscopic, Short-wire				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
BKM0604CW	6	4	8.5	180
BKM0605CW		5		
BKM0606CW		6		
BKM0607CW		7		
BKM0608CW		8		
BKM0609CW		9		
BKM0610CW		10		
BKM0612CW		12		
BKM0804CW		8		
BKM0805CW	5			
BKM0806CW	6			
BKM0807CW	7			
BKM0808CW	8			
BKM0809CW	9			
BKM0810CW	10			
BKM0812CW	12			
BKM1004CW	10		4	8.5
BKM1005CW		5		
BKM1006CW		6		
BKM1007CW		7		
BKM1008CW		8		
BKM1009CW		9		
BKM1010CW		10		
BKM1012CW		12		

Endoscopic				
Bumpy with a long Suture				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
BK0604CW2	6	4	8.5	180
BK0606CW2		6		
BK0608CW2		8		
BK0804CW2	8	4	8.5	180
BK0806CW2		6		
BK0808CW2		8		
BK1004CW2	10	4	8.5	180
BK1006CW2		6		
BK1008CW2		8		