


# Niti-S Enteral Colonic Uncovered Stent [D-Type]

for intrinsic and/or extrinsic malignant colonic stricture



## Features

- **Unfixed cell with weaving construction**
  - Conformability to prevent migration
  - Combination of radial and axial force to maintain patency in tortuous anatomy
- **Radiopaque marker: 4 (four) at both ends & 2 (two) in the middle**

# Ordering Information

TTS (Through The Scope)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
CDT1806	18	6	10	220
CDT1808		8		
CDT1810		10		
CDT1812		12		
CDT1814		14		
CDT1815		15		
CDT2006	20	6	10	220
CDT2008		8		
CDT2010		10		
CDT2012		12		
CDT2014		14		
CDT2015		15		
CDT2206	22	6	10	220
CDT2208		8		
CDT2210		10		
CDT2212		12		
CDT2214		14		
CDT2215		15		
CDT2406	24	6	10	220
CDT2408		8		
CDT2410		10		
CDT2412		12		
CDT2414		14		
CDT2415		15		
CDT2606	26	6	10.5	220
CDT2608		8		
CDT2610		10		
CDT2612		12		
CDT2614		14		
CDT2615		15		
CDT2806	28	6	10.5	220
CDT2808		8		
CDT2810		10		
CDT2812		12		
CDT2814		14		
CDT2815		15		

TTS (Through The Scope)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
CDT1806-12	18	6	10	220
CDT1808-12		8		
CDT1810-12		10		
CDT1812-12		12		
CDT1814-12		14		
CDT1815-12		15		
CDT2006-12	20	6	10	220
CDT2008-12		8		
CDT2010-12		10		
CDT2012-12		12		
CDT2014-12		14		
CDT2015-12		15		
CDT2206-12	22	6	10	220
CDT2208-12		8		
CDT2210-12		10		
CDT2212-12		12		
CDT2214-12		14		
CDT2215-12		15		
CDT2406-12	24	6	10	220
CDT2408-12		8		
CDT2410-12		10		
CDT2412-12		12		
CDT2414-12		14		
CDT2415-12		15		
CDT2606-12	26	6	10	220
CDT2608-12		8		
CDT2610-12		10		
CDT2612-12		12		
CDT2614-12		14		
CDT2615-12		15		
CDT2806-12	28	6	10	220
CDT2808-12		8		
CDT2810-12		10		
CDT2812-12		12		
CDT2814-12		14		
CDT2815-12		15		

# Ordering Information

OTW (Over The Wire)				
Distal Release				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
CD1806	18	6	16	70
CD1808		8		
CD1810		10		
CD1812		12		
CD1814		14		
CD1815		15		
CD2006	20	6	16	70
CD2008		8		
CD2010		10		
CD2012		12		
CD2014		14		
CD2015		15		
CD2206	22	6	16	70
CD2208		8		
CD2210		10		
CD2212		12		
CD2214		14		
CD2215		15		
CD2406	24	6	16	70
CD2408		8		
CD2410		10		
CD2412		12		
CD2414		14		
CD2415		15		
CD2606	26	6	18	70
CD2608		8		
CD2610		10		
CD2612		12		
CD2614		14		
CD2615		15		
CD2806	28	6	18	70
CD2808		8		
CD2810		10		
CD2812		12		
CD2814		14		
CD2815		15		
CD3006	30	6	18	70
CD3008		8		
CD3010		10		
CD3012		12		
CD3014		14		
CD3015		15		

OTW (Over The Wire)				
Proximal Release				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
CDP1806	18	6	16	70
CDP1808		8		
CDP1810		10		
CDP1812		12		
CDP2006	20	6	16	70
CDP2008		8		
CDP2010		10		
CDP2012	22	12	16	70
CDP2206		6		
CDP2208		8		
CDP2210		10		
CDP2212	24	12	16	70
CDP2406		6		
CDP2408		8		
CDP2410		10		
CDP2412	26	12	18	70
CDP2606		6		
CDP2608		8		
CDP2610		10		
CDP2612	28	12	18	70
CDP2806		6		
CDP2808		8		
CDP2810		10		
CDP2812	30	12	18	70